
Hydrogen sulphide in liquids and gases

Conform to standard DIN 38405-27







sulfimax GX

Hydrogen sulphide in liquids and gases

Conform to standard DIN 38405-27

Product description

The **Sulfimax GX Go** determines hydrogen sulphide and volatile sulphides in liquids and gases.

The sample can be dispensed directly by suringe without pretreatment.

By effective gas extraction the H₂S is completely expelled from the sample. Interferences due to the sample matrix practically do not occur.


The released H₂S gas is conducted to the highly sensitive sensor, which detects H₂S in the range of 0.01 to 10,000 ppm. A typical measurement takes 5 min, depending on the sample composition.

If the **Sulfimax GX Go** is extended with the optionally available H₂S Headspace Module, solid and pasty samples can also be measured.

The Sulfimax GX Lab fulfils the requirements of the standard DIN 38405-27: German standard methods for the examination of water, waste water and sludge - Anions (group D) - Part 27:

> Determination of sulphide by gas extraction method (D 27).





Sulfimax GX Go - compact version with small footprint

Applications

- Water, drinking water, surface water
- Municipal wastewater
- Industrial wastewater
- Monitoring of landfill-leachate
- Gas analysis (e. g. LNG, LPG)
- H₂S in hydrocarbon mixtures
- Investigation of technical and pharmaceutical products (e.g. storage stability)
- Quality management

Advantages


- Analysis of the original sample
- No sample preparation
- · Minimized cross sensitivity through the indirect method
- Suitable for on-site use




Sulfimax GX Go - interior view

Features and Results

- Complete separation of H₂S from the sample by effective gas extraction
- Simple calibration
- Software: simple, clear, intuitive
- Dosing manually by syringe
- Definition of own methods for device control
- Typical measuring time 5 min (depending on sample properties)



Functional scheme of the analysis system



Typical measurement - automatic peak analysis / interpretation




Table of results of a multi measurement

Technical specifications

Measuring range: 0.01 ... 10,000 ppm (dependent on sample volume)

Resolution: 0.1 µg abs., output signal linear

Measuring duration: 2 ... 10 min (dependent on the sample), usually 5 min

Sample volume: $0.01 \dots 20 \text{ mL}$ Gas flow: Up to 50 L/h

Power supply: 230 V/50 Hz, 115 V/60 Hz

Power input: 150 W

Dimensions: $260 \times 310 \times 300 \text{ mm} (W \times D \times H)$

Weight: 8 kg

Data connection: RS 232 / USB (with converter)

Device control: PC software (PC not included in the scope of delivery)



Laboratory version **Sulfimax GX Lab** with autosampler for liquids

H₂S Headspace Module

Extension module for solid and pasty samples

The determination of volatile hydrogen sulphide (H₂S) in solid and pasty samples is easily possible with this module.

It is connected directly to the selective **Sulfimax GX Go** and can be operated by anyone.

The sample is heated isothermally without pretreatment in a sealed headspace vessel and analyzed.

Applications:

- Solid samples, e. g. elemental sulfur, sludge, bitumen
- Liquid samples like waste water with sludge particles
- Pasty samples
- Soil samples and waste



Sulfimax GX Go with connected H₂S Headspace Module

Book your online demo in the ECH Studio

ECH Scientific have a state of the art laboratory fitted with online presentation capabilities, allowing us to bring product demonstrations live and in full HD, with multiple camera angles and software sharing capabilities enabling us to deliver a full demo experience remotely. Please contact info@echscientific.com to book your session.

ECH Elektrochemie Halle GmbH

Otto-Eißfeldt-Str. 8 D-06120 Halle (Saale)

Germany

Tel.: +49 (0) 345 279570-0 Fax: +49 (0) 345 279570-99 **ECH Scientific Limited**

Building 69, Wrest Park, Silsoe Bedfordshire, MK45 4HS

United Kingdom

Tel.: **+44 (0) 1525 404747** Fax: **+44** (0) 1525 404848

Email: info@echscientific.com • www.ech.de • www.aquamaxkf.com



the ECH advantage