



# Specification Sheet

---

<b>Product Name</b>	PENTANE
<b>Item Number</b>	7801-7-40

---

<b>Test Description</b>	<b>Limit</b>
Assay (by GC)	Min. 99.0 %
<b>Maximum allowable</b>	
Water (coulometric KF)	0.01 %
Residue after evaporation	3 ppm
<b>UV Abs (1.00cm path vs. water)</b>	
UV @ 190nm	1.00 AU
UV @ 200nm	0.40 AU
UV @ 220nm	0.06 AU
UV @ 245nm	0.01 AU
UV @ 254nm	0.01 AU
UV @ 300nm	0.01 AU
UV @ 350nm	0.01 AU
UV @ 400nm	0.01 AU
UV Cutoff	190 nm

---

Caledon Laboratories, Ltd. believes that the information contained herein is reliable and accurate. Caledon makes no warranty thereto, and expressly disclaims all liability for reliance thereon. Such information is solely for your consideration.

40 Armstrong Avenue, Georgetown, Ontario L7G 4R9 Phone (905) 877-0101, Fax (905) 877-6666