




INTER DYNE  
SYSTEMS

Inter Dyne Systems, Inc.  
676 Ellis Road  
Norton Shores, MI 49441  
Ph: 877.4ss.labs  
www.interdynesystems.com

## Pegboard Accessories

## Cylinder Yolk

Model# CY-12



### Features

- 14"w x .500"d
- Mounts on pegboard using 3 SP-4 support pegs
- Autoclavable
- Constructed from 16 gauge, type 304 stainless steel, polished to a #4 finish
- Shipped with 3 SP-4 pegs & 2 HP-1 hole plugs to mount on pegboard

### Product Dimensions

