




INTER DYNE
SYSTEMS

Inter Dyne Systems, Inc.
676 Ellis Road
Norton Shores, MI 49441
Ph: 877.4ss.labs
www.interdynesystems.com

Pegboard Accessories

Funnel Rack

Model# FR-12



Features

- 14"w x .375"h x 5.07"d
- Mounts on pegboard using 3 SP-4 support pegs
- Autoclavable
- Constructed from 20 gauge, type 304 stainless steel, polished to a #4 finish
- Shipped with 3 SP-4 pegs & 3 HP-1 hole plugs to mount on pegboard

Product Dimensions

