

Specification Sheet

Product Name TOLUENE

Item Number 9201-7-40

Test Description	Limit
------------------	-------

Also meets ACS Reagent specifications

Assay (by GC) Min. 99.8 %

Maximum allowable

Colour 10 APHA
Substances darkened by sulphuric acid Passes Test
Sulphur compounds (as S) 0.003 %
Residue after evaporation 3 ppm
Water (coulometric KF) 0.02 %

UV Abs (1.00cm path vs. water)

UV @ 285nm 1.00 AU
UV @ 288nm 0.40 AU
UV @ 293nm 0.20 AU
UV @ 300nm 0.10 AU
UV @ 310nm 0.05 AU
UV @ 330nm 0.02 AU
UV @ 335nm 0.02 AU
UV @ 350nm 0.01 AU
UV @ 400nm 0.01 AU
UV Cutoff 285 nm

Caledon Laboratories, Ltd. believes that the information contained herein is reliable and accurate. Caledon makes no warranty thereto, and expressly disclaims all liability for reliance thereon. Such information is solely for your consideration.

40 Armstrong Avenue, Georgetown, Ontario L7G 4R9 Phone (905) 877-0101, Fax (905) 877-6666