



# Specification Sheet

---

<b>Product Name</b>	TOLUENE
<b>Item Number</b>	9200-1-40

---

<b>Test Description</b>	<b>Limit</b>
<b>Meets ACS Specifications</b>	
Assay (by GC, corrected for water)	Min. 99.5 %
<b>Maximum allowable</b>	
Colour	10 APHA
Residue after evaporation	0.001 %
Substances darkened by sulphuric acid	Passes Test
Sulphur compounds (as S)	0.003 %
Water (coulometric KF)	0.03 %

---

Caledon Laboratories, Ltd. believes that the information contained herein is reliable and accurate. Caledon makes no warranty thereto, and expressly disclaims all liability for reliance thereon. Such information is solely for your consideration.

40 Armstrong Avenue, Georgetown, Ontario L7G 4R9 Phone (905) 877-0101, Fax (905) 877-6666