



Specification Sheet

Product Name	XYLENES
Item Number	9800-1-40

Test Description	Limit
Meets ACS Specifications	
Assay (Xylene isomers and EB corrected for water)	Min. 98.5 %
Assay (Ethyl benzene)	Max. 25 %
Maximum allowable	
Colour	10 APHA
Residue after evaporation	0.002 %
Substances darkened by sulphuric acid	Passes Test
Sulphur compounds (as S)	0.003 %
Water (coulometric KF)	0.05 %

Caledon Laboratories, Ltd. believes that the information contained herein is reliable and accurate. Caledon makes no warranty thereto, and expressly disclaims all liability for reliance thereon. Such information is solely for your consideration.

40 Armstrong Avenue, Georgetown, Ontario L7G 4R9 Phone (905) 877-0101, Fax (905) 877-6666