

Eco smartFunnel™

Paper Funnel

A clean transfer every time

- Supplied flat-pack to minimize storage
- With a minimal amount of finger pressure the funnel takes its working shape
- Maximum flexibility – cut the stem to any size or shape to meet with your needs
- Fits a wide variety of openings and can be used with virtually any liquid* or powder



* Not recommended for hot liquids

Pack of 25

Item No.	Capacity mL	Top Dia. in ²	Top Dia. cm ²	Bottom Dia. in ²	Bottom Dia. cm ²	Height in	Height cm
120219	186	2.4	6	0.3	0.6	5	12.8

Eco smartFunnel™ is a trademark of LevGo Inc.

Lead Rings

Stabilize vessels and labware on the bench or in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit



Item No.	Shape	Color	Flask Sizes mL	ID in	ID cm	OD in	OD cm	Weight g
HS8881A	C	Blue	125–500	1.9	4.2	3	7.5	215
HS8882A	O			1.9	4.8	3	7.7	210
HS8881B	C	Blue	250–1000	2	5.1	3.5	9	465
HS8882B	O			2	5.1	3.5	9	490
HS8881C	C	Blue	500–2000	2.1	5.4	3.9	9.8	595
HS8882C	O			2.2	5.7	3.9	9.8	610
HS8881D	C	Blue	1000–4000	2.6	6.6	4.4	11.3	890
HS8882D	O			2.8	7	4.5	11.5	935

*Weights may vary by +/- 3%

Safety Bottle Carriers

Secure and convenient bottle transport

- Ergonomic handle prevents hand fatigue and related accidents when carrying heavier loads
- Handle swings down to either side for easier loading and unloading
- Centering ring stabilizes bottle in the center of carrier for safer, steadier transport
- Carriers come in four colors with designated labeling areas for easy identification
- Made of chemical resistant thermoplastic elastomer - PP blend



Available in 4 Colors
for easy identification and organization

Safety Bottle Carrier, 4.5 L

Item No.	Color	Dimensions (Dia. x H) in	Dimensions (Dia. x H) cm
120911	Red	7.2 x 10.5	18.3 x 26.7
120912	Black		
120913	Blue	(16 in height, handle extended)	(40.6 cm height, handle extended)
120914	Green		

Safety Bottle Carrier, 1.5 L

Item No.	Color	Dimensions (Dia. x H) in	Dimensions (Dia. x H) cm
120915	Red	4.5 x 6.5	11.4 x 16.5
120916	Black		
120917	Blue	(10 in height, handle extended)	(25.4 cm height, handle extended)
120918	Green		



Movable Handle
swings down to allow
clear access to vessel

