



CENTRIFUGE TUBE OVERVIEW




Demand a Higher Standard

Tel: (713) 926-2623
Fax: (713) 926-7736


*Petroleum Testing Equipment
Suppliers Worldwide Since 1930*

Centrifuge Tube Types:



| Finger (A) | Pear-Shaped (B) | Long Cone (C) | Short Cone (D) | Capillary (E) |
|---|---|---|---|--|
| <ul style="list-style-type: none"> Used with hand-crank or Mini-Melton API 2542 | <ul style="list-style-type: none"> Extended capillary tip for accurate readings API 2542, ASTM D96 (API 10.4) | <ul style="list-style-type: none"> Laboratory or refinery (Benchmark models) ASTM D4007 (API 10.3), ASTM D96 (API 10.4) | <ul style="list-style-type: none"> Transport models with limited space ASTM D4007 (API 10.3), ASTM D96 (API 10.4) | <ul style="list-style-type: none"> Best for refined fuel products ASTM D2709 |

Available Markings, Tops & Threads:



Standard (100-ml)

Detailed (100-ml)

Standard (200%)

Detailed (200%)

Cork

Wide-Neck Twist Cap

Rubber Stopper

Standard Twist Cap

Cork Toolled Top

Wide-Neck Threaded

Rubber Toolled Top

Standard Threaded