

RUGGED, ECONOMICAL WIPING SEAL DISPENSERS

seripettor® & seripettor® pro

The seripettor® and seripettor® pro bottle-top dispensers feature an innovative replaceable piston cylinder cartridge that allows easy renewal of wearing parts. They are an excellent choice for labs dispensing crystallizing media, or labs where maintenance responsibility is shared. They also are an economical alternative

for customers who do not need the safety features or variety of accessories available for the Dispensette® S.

A spring-loaded, semi-automatic design greatly facilitates repetitive sample preparation.

seripettor®

The seripettor® is suitable for many common liquids found in Biological laboratories: bases; buffers; culture media; hydrogen peroxide; or dilute, non-oxidizing acids. Even sodium azide solutions can be dispensed, something that destroys more expensive dispensers.

seripettor® pro

Choose the seripettor® pro for more aggressive media, like DMSO, acetone, or dilute hydrochloric or nitric acids. An amber cylinder sleeve and opaque outer housing helps protect light-sensitive media.

Each seripettor® or seripettor® pro bottle-top dispenser includes:

- Spare dispensing cartridge
- Discharge tube
- Filling tube
- Three PP adapters (GL 45/33, GL 45/38 & GL 45/S40)
- Operating manual
- One-year warranty

Description	BrandTech® Cat. No.
seripettor®	
0.2-2 mL	4720120
1-10 mL	4720140
2.5-25 mL	4720150
seripettor® pro	
0.2-2 mL	4720420
1-10 mL	4720440
2.5-25 mL	4720450

Product appearance, catalog numbers, prices, specifications, and technical information are subject to change without notice.



For assistance in selecting the correct dispenser, ask your BrandTech® Sales Representative, or visit www.brandtech.com.