

Areas of application / Chemical selection list

■ seripettor® dispenser

■ seripettor® pro dispenser

Reagent	seripettor®	seripettor® pro®	Reagent	seripettor®	seripettor® pro®
Acetaldehyde		+	Glycol (Ethylene glycol)	+	+
Acetic acid, 5%	+	+	Glycolic acid, 50%	+	+
Acetic acid, 96%		+	Hexanoic acid	+	+
Acetic acid (glacial), 100%		+	Hexanol		+
Acetone		+	Hydriodic acid	+	+
Acetonitrile		+	Hydrobromic acid		+
Acetophenone	+		Hydrochloric acid, 37%		+
Acetylacetone	+	+	Hydrogen peroxide, 35%	+	
Acrylic acid	+	+	Isoamyl alcohol		+
Acrylonitrile		+	Isobutanol	+	+
Adipic acid	+	+	Isopropanol (2-Propanol)	+	+
Agar (60°C)	+		Lactic acid	+	+
Allyl alcohol	+	+	Methanol	+	+
Aluminium chloride	+	+	Methyl benzoate		+
Amino acids	+	+	Methyl ethyl ketone		+
Ammonia 30%	+	+	Methyl propyl ketone		+
Ammonium chloride	+	+	Mineral oil (Engine oil)		+
Ammonium fluoride	+	+	Monochloroacetic acid		+
Ammonium sulfate	+	+	Nitric acid, 10%		+
Amyl alcohol (Pentanol)	+	+	Oxalic acid	+	+
n-Amyl acetate		+	Perchloric acid		+
Aniline		+	Phenol		+
Barium chloride	+	+	Phosphoric acid, 85%		+
Benzaldehyde		+	Piperidine		+
Benzyl alcohol		+	Potassium chloride	+	+
Benzylamine		+	Potassium dichromate	+	+
Benzylchloride		+	Potassium hydroxide	+	+
Boric acid, 10%	+	+	Potassium hydroxide in ethanol	+	+
Butanediol	+	+	Potassium permanganate	+	+
1-Butanol		+	Propionic acid	+	+
Butylamine		+	Propylene glycol (Propanediol)	+	+
n-Butyl acetate		+	Pyridine		+
Calcium carbonate	+	+	Pyruvic acid	+	+
Calcium chloride	+	+	Salicylaldehyde		+
Calcium hydroxide	+	+	Salicylic acid	+	+
Calcium hypochlorite		+	Silver acetate	+	+
Chloroacetaldehyde, 45%		+	Silver nitrate	+	+
Chloroacetic acid		+	Sodium acetate	+	+
Chromic acid, 50%		+	Sodium chloride	+	+
Copper sulfate	+	+	Sodium dichromate	+	+
Cumene (Isopropyl benzene)		+	Sodium fluoride	+	+
Diethylene glycol	+	+	Sodium hydroxide, 30%	+	+
Dimethyl sulfoxide (DMSO)		+	Sodium hypochlorite (active Chlorine approx. 10%)		+
Dimethylaniline		+	Sulfuric acid, 10%	+	+
Ethanol	+	+	Tartaric acid		+
Formaldehyde, 40%	+	+	Urea	+	+
Formamide	+	+	Zinc chloride, 10%	+	+
Formic acid, 100%		+	Zinc sulfate, 10%	+	+
Glycerol	+	+			

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BrandTech. Status as of: 0114/9.

For the most current version of this chart, visit our website www.brandtech.com.

NOTE: The seripettor® & seripettor® pro bottle top dispensers cannot be used with hydrofluoric acid (HF).